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01	MODELS.....	36
01.1	Products.....	36
01.2	Accessories.....	36
02	INTRODUCCIÓN.....	37
02.1	Using this guide.....	37
02.2	Legend of symbols.....	37
02.3	Service request.....	38
02.4	Life span.....	38
02.5	Demolition.....	39
02.6	Reference standards.....	39
02.7	Labelling.....	40
03	INTENDED USE.....	41
03.1	Intended use.....	41
03.2	Contraindications and side effects.....	41
03.3	Physical requirements of the operators.....	41
04	WARNINGS.....	42
04.1	General warnings.....	42
04.2	Specific warnings.....	44
04.3	Residual risk.....	45
05	DESCRIPTION OF THE PRODUCT.....	46
05.1	Main components.....	46
05.2	Technical data sheet.....	48
05.3	Main components (accessories).....	49
05.4	Technical data sheet (accessories).....	50
06	OPERATION.....	51
06.1	Transport and storage.....	51
06.2	Preparation.....	51
06.3	Functioning.....	52
07	GENERAL MAINTENANCE.....	56
07.1	Cleaning.....	56
07.2	Lubrication.....	56
07.3	Wear areas.....	57
07.4	Maintenance program.....	57
08	SPARE PARTS.....	59
09	TRAINING REGISTER.....	61
10	MAINTENANCE REGISTER.....	62
11	LEGAL NOTICES.....	63
12	PRODUCT WARRANTY.....	64

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REV. 2023/10

01.1. Products

REFERENCE	UDI CODE	TRADE NAME
PS1900-00000	8437001519149	EVACUATION CHAIR WITH TRACKS PS-190
PS1910-00000	8437001519156	EVACUATION CHAIR WITH TRACKS PS-191
PS1950-00000	8436601870742	EVACUATION CHAIR WITH TRACKS PS-195 (RED)
PS1960-00000	8436601870643	EVACUATION CHAIR WITH TRACKS PS-196 (WHITE)
PS1970-00000	8436601870759	EVACUATION CHAIR LARGE WHEELS WITH TRACKS (BLUE)
PS1980-00000	8436601870728	EVACUATION CHAIR LARGE WHEELS WITH TRACKS (ORANGE)

01.2. Accessories

REFERENCE	UDI CODE	TRADE NAME
PA5780-00000	8437001519217	LOCK SET PA-578
PA5790-00000	8436601870896	LOCK SET PA-579

02 INTRODUCTION

REV. 2023/10

02.1. Using this guide

The manual provides instructions for use and maintenance of the product, as well as technical aspects, operation, spare parts and safety.

It is recommended before the operation of the product to read carefully this manual in order to avoid damages caused by a misuse, regardless of the level of experience acquired in the past with similar devices.

Do not lose this document. It should be accessible to any doubt that could appear by medical personnel. In case the device has some instructions for use of a different product than the one received, it is necessary to contact the manufacturer immediately before using the product.

Remember that a good use and maintenance are necessary for the proper operation of the product.

Each product includes an identifying adhesive label with the serial number and model. Keep these numbers so that you can refer them to the dealer if necessary.

02.2. Legend of symbols

SYMBOL	EXPLANATION / DESCRIPTION
	MANUFACTURER symbol. This symbol is accompanied by the name and address of the manufacturer, adjacent to the symbol (PRODUCTOS METÁLICOS DEL BAGÉS S.L., Ctra. C-16 Km 59.5, 08650 Sallent (Barcelona)).
	Indicates the manufacturer's reference number to identify the medical device. PROMEBA, S.L. uses this symbol to set each internal reference for each configuration and business variant
	Indicates the manufacturer's serial number to identify a specific medical device.
	Indicates the manufacturing date. The symbol must be accompanied by a manufacturing date (yyyy-mm), adjacent to the symbol.
	It is placed to inform that the product is a "Medical Device"
	CE symbol without the intervention of a Notified Organism, as a Medical Device classified as Class I according to the EU Regulation 2017/745 on Medical Devices.
	Symbol for the Unique Device Identifier.
	Symbol "See instructions for use or operating instructions"
	Symbol "caution". This symbol is placed to warn of the need for the user to refer to important precautionary information in the operating instructions, such as warnings and cautions not otherwise found on the label.
	Waste symbol for electrical and electronic equipment "WEEE" (Waste Electrical and Electronic Equipment). Recycle: Electronic equipment. DO NOT THROW GARBAGE
	"Dangerous voltage" symbol. To warn of electricity

02 INTRODUCTION

REV. 2023/10


Symbol "Caution". For a general warning.


Warning, crushing of hands

"Double insulation" symbol. The protection of appliances marked with this symbol is ensured by double insulation and does not require a safety electrical connection to earth (ground).

Cot duty cycle: 16.7% (less than 60 seconds on, more than 300 seconds off)

Indicates the temperature limits. The upper and lower temperature limits should be indicated adjacent to the horizontal lines.

02.3. Service request

For information of the correct interpretation of the instruction manual, the use, maintenance, installation and restoration of the product, please contact Promeba customer service: T. 93 837 12 00, email: promeba@promeba.com or write to PROMEBA, S.L. - Ctra C-16 Km 59.5 - 08650 Sallent (Barcelona) - SPAIN.

Whenever you want to communicate any matter related to a product, you will be asked for the product reference (**REF**) and the serial number (**SN**). It is strongly recommended to keep the label containing this information and found on the original packaging of the product.

02.4. Life span

If used as described in the following instructions, this device has a useful life of 10 years from the date of purchase.

This useful life can be extended with annual reviews carried out by the manufacturer, which uses specialized and authorized internal and external technicians.

In case these annual checks are not carried out, the device must be disposed of according to the information in paragraph 02.5 and the manufacturer must be notified.

Only the manufacturer or an authorized center can extend the life of the device, if it meets the safety and legal requirements.

Promeba, S.L. will not accept any responsibility for malfunction or damage caused by the use of devices that have not been checked by the manufacturer or authorized center, or that have exceeded the maximum allowed useful life.

02 INTRODUCTION

REV. 2023/10

02.5. Demolition

When the devices are no longer suitable for use, if they have not been contaminated by any particular agent, they can be disposed of as electric waste, otherwise, follow the current demolition regulations.

WARNINGS FOR THE CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/EU WEEE:

At the end of his life, the product must not be disposed as household waste. Can be taken to special recycling centers provided by local government, or return it to the dealer on purchase of a new device of the same type and used for the same functions. Dispose of the product separately avoids possible negative consequences for the environment and human health resulting from inappropriate disposal and allows to recover the materials in order to obtain significant savings in energy and resources. The symbol on the label indicates separate collection of electrical and electronic equipment.

02.6. Reference standards

REFERENCE	TITLE OF DOCUMENT
EU Regulation 2017/745	EU Regulation on Medical Devices
EN 1789:2007+A2:2014	Medical vehicles and their equipment - Road ambulances
UNE-EN 1865-4:2012	Patient handling equipment used in road ambulances - Part 4: Foldable patient transfer chair
EN 60601-1:2006/A1:2013 IEC 60601-1:2005/A1:2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
EN 61000-6-1:2007 IEC 61000-6-1:2005	Electromagnetic compatibility (EMC) - Part 6-1: Generic standards - Immunity for residential, commercial and light-industrial environments
EN 61000-6-3:2007+A1:2012 IEC 61000-6-3:2006+A1:2010	Electromagnetic compatibility (EMC) -- Part 6-3: Generic standards - Emission standard for residential, commercial and light-industrial environments


As a distributor or end user of the products manufactured and/or distributed by Promebea, S.L., you are strictly obliged to know the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including the regulations relating to technical specifications and/or safety requirements) and, therefore, understand the necessary requirements to ensure compliance of the products themselves with all the legal requirements of the territory.

02 INTRODUCTION

REV. 2023/10

02.7. Labelling

Each product incorporates an identification label, placed on the device itself and/or on the original packaging of the product. It must never be removed or covered. This label includes the data of the company which produced the device (🏭), the product code / reference (REF) of the product, the manufacturing date (📅), the serial number (SN), the UDI-PI code (UDI) and the CE distinctive (CE). Please keep these numbers so you can inform the dealer if necessary.

The UDI-PI code is made up of the following fields:

- (01)08437001519149: UDI code of the product.
- (11)231009: it gives information about the batch release date of the product following the format *aammdd*, where *aa* is the batch release year, *mm* is the month of the batch release and *dd* is the day of the batch release.
- (21)01564123001: identifies the serial number of the product giving information about the manufacturing order (0156) including the week when the product was produced (4123) and a correlative number depending on the units that has been produced in that production order (001).

03 INTENDED USE

REV. 2023/10

03.1. Intended use

The product ELECTRIC CHAIR is a device indicated to effortlessly transport patients (over 40kg) in hospital and pre-hospital environments thanks to its electric transport function when going up and down stairs.

03.2. Contraindications and side effects

Any contraindication or side effect have been resulted from the clinical analysis of the ELECTRIC CHAIR family of products.

03.3. Physical requirements of the operators

Promeba electric chairs are intended exclusively for professional use. Each operator must be trained to transport patients safely and efficiently. Do not allow untrained persons to assist during use of the product, as they may cause injury to others or themselves. The device is intended to be used by medical transport technical personnel or ambulance technicians who have knowledge and experience with patient transport.

Operators must be able to provide the necessary assistance to the patient.

Users must be able to safely lift and manipulate the weight of the assembly consisting of the chair, the patient, and any other equipment used with the device.

For patient loading techniques, for particularly heavy patients, for interventions on steep terrain or in particular and unusual circumstances, the presence of several operators is recommended (not just 1 as expected in standard conditions).

Operators must be able to evaluate any anomaly of the system/assembly, communicating the problem to those responsible and interrupting the use of the devices if necessary to ensure patient safety.

The capacities of the various operators must be considered before determining his role in the employment of the board.

04 WARNINGS

REV. 2023/10

04.1. General warnings

1. The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
2. Establish a maintenance and periodic control program and identify a reference manager. The person in charge of the ordinary maintenance of the device must guarantee the basic requirements provided by the manufacturer in these instructions for use. It is recommended that the subject check at least every 6 months if there are updated user manuals available or if there are any changes related to his / her product. This information is freely available on the website <http://promeba.com/>.
3. Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register, which will certify the eligibility of the operators to use the Promeba, S.L. device, has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
4. Promeba, S.L. is always at your disposal to plan trainings on products.
5. Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
6. If the instructions belong to another device and not the device received, inform the manufacturer immediately and avoid use of the device.
7. In the case of any doubts as to the correct interpretation of the instructions, please contact Promeba, S.L. for any necessary clarifications.
8. Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
9. Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and/or of the user are detected, the device must be immediately removed from service and the manufacturer must be contacted.
10. If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
11. Use of the device in anyway other than described in this manual is forbidden.
12. Do not alter or modify in any way the device; any such interference could cause malfunctions and injury to the patient and/or rescuer.
13. The device must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the device itself; moreover CE certification and product warranty will be considered void.

04 WARNINGS

REV. 2023/10

14. Use only original or PROMEBA approved components/spare parts and/or accessories to carry out all operations without causing alterations or modifications to the device; Otherwise, any responsibility for the operation is declined.
15. Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
16. Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users and any other comments.
17. When the device is being used, the assistance of qualified staff must be guaranteed.
18. The product must not be exposed or come into contact with thermal sources of combustion and flammable agents, but must be stored in a dry, cool place, protected from light and sunlight.
19. Store and transport device in its original packaging.
20. Do not store the product under other more or less heavy materials, which can damage the structure of the product.
21. Attention: laboratory testing, post production tests and instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
22. As a distributor or end users of products manufactured and/or marketed by Promeba, S.L., you are strictly required to have a basic knowledge of any legal requirements existing in the goods final destination country applying to the devices contained in this supply (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
23. Promptly notify PROMEBA regarding any revisions to be made by manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
24. Act with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant user manual.
25. Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary actions can be promptly taken.
26. Be aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore we expressly disclaim any responsibility and/or liability for your non-compliance with the present "Regulatory provisions".
27. Follow internal procedures and protocols approved by your organization.

04 WARNINGS

REV. 2023/10

28. Follow approved Emergency Medical Service procedures for patient immobilization and transportation.

04.2. Specific warnings

1. Use only accessories/spare parts that are original or approved by Promeba, S.L. when performing any operation so as not to cause any alteration or modification to the device, otherwise we do not assume any responsibility for the correct operation or resulting damage of the device to the patient or the operator and the guarantee will be considered void, as well as compliance with the Regulation 2017/745.
2. Always respect the maximum capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the load of the total weight on the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
3. Inside the ambulance, place and adjust the device taking care not to obstruct other parts or other equipment. It is recommended not to open the equipment inside the ambulance, but rather to do it outside the vehicle.
4. Never leave the patient unassisted on the device, because he may be injured.
5. Do not use bleach to disinfect the device. Use a hydroalcoholic based disinfectant and wash with water.
6. The device and all its components, after washing, should be allowed to dry completely before storing.
7. Avoid contact with sharp objects.
8. Do not use the device if it is pierced, deformed or shows signs of loss of integrity.
9. Do not operate in case the weight has not been distributed correctly.
10. Condensation, water, ice and dust accumulations can affect the correct operation of the device, making it unpredictable.
11. Always immobilize the patient at least with the belts supplied by the manufacturer during any movement of the stretcher, including loading and unloading to the ambulance.
12. Do not use the device as a support for climbing or leaning on and reaching high points in the ambulance. The product may be damaged and the technician or patient may be harmed.
13. Disinfection activities must be carried out in accordance with the validated cycle parameters, indicated in the specific technical standards.
14. Promeba S.L. shall not be liable if unauthorized / incompatible stretcher fastener is used.
15. Do not leave the batteries uncharged for more than 2 MONTHS. They could break down.

04 WARNINGS

REV. 2023/10

16. Incorrect loading of the chair can cause damage and/or injuries. Use only batteries approved and authorized by Promeba S.L. and recharge only with chargers approved and authorized by Promeba S.L.
17. Do not use the chair with the track carriage without being properly locked as this may cause damage and/or injury. Make sure you have properly locked the track carriage before moving the patient.
18. Do not exceed the maximum capacity allowed by the chair. Exceeding it may cause deformations to the product or harm to the patient.
19. Do not lubricate the chair when it is connected to the power.
20. In repeated up and down maneuvers of up to 10 uninterrupted cycles, a mandatory stop must be made for at least 15 seconds so that the electric chair's motor-reduction system recovers.
21. Always use both hands to transport the chair to avoid the risk of tipping over.
22. Pay full attention to the team and the patient during the performance of the service.

04.3. Residual risk

The residual risks listed below have been identified only with reference to the intended use of the device:

1. A patient without assistance is at risk of injury. Never leave the patient alone and always help them while using the chair.
2. The operation carried out by unqualified personnel could cause harm to both the patient and the personnel performing the operation.

05.1. Main components

PS1900-00000 EVACUATION CHAIR WITH TRACKS PS-190

PS1910-00000 EVACUATION CHAIR WITH TRACKS PS-191

N°	DESCRIPTION OF THE COMPONENTS	N°	DESCRIPTION OF THE COMPONENTS
1	Armrests	7	Rear wheel with brake
2	Belts	8	Front swivel wheels
3	Telescopic handles	9	Control panel
4	Footrest	10	Folding system with safe anti-fold
5	Telescopic backrest	11	Track system with anti-slide
6	Foldable transport handles		

PS1950-00000 EVACUATION CHAIR WITH TRACKS PS-195 (RED)
PS1960-00000 EVACUATION CHAIR WITH TRACKS PS-196 (WHITE)

Nº	DESCRIPTION OF THE COMPONENTS	Nº	DESCRIPTION OF THE COMPONENTS
1	Armrests	7	Rear wheel with brake
2	Belts	8	Front swivel wheels
3	Telescopic handles	9	Control panel
4	Footrest	10	Folding system with safe anti-fold
5	Telescopic backrest	11	Track system with anti-slide
6	Foldable transport handles	12	Seat and backrest in thermoformed ABS

05.2. Technical specifications / data sheet

Models	PS-190	PS-191	PS-195	PS-196
Length	2090mm	2060mm	685/1180mm	685/1180mm
Width	550mm	550mm	550mm	550mm
Height	170mm	70mm	1090mm	1090mm
Weight	30kg	30kg	34kg	34kg
Max. Load	180kg	180kg	180kg	180kg
Color	Rojo	Blanco	Rojo	Blanco
Ø front wheel	75mm	75mm	150mm	150mm
Ø rear wheel	125mm	125mm	200mm	200mm
Batery duration	1.5h	1.5h	1.5h	1.5h

05.3. Main parts (accessories)

PA5780-00000 LOCK SET PA-578	
PA5790-00000 LOCK SET PA-579	
N°	DESCRIPTION OF THE COMPONENTS
1	Belt fastening system
2	Centering platform
3	Wheels slots

05.4. Technical specifications / data sheet (accessories)

Modelss	PA5780-00000	PA5790-00000
Length	550mm	550mm
Width	219mm	219mm
Height	44 / 218mm	44 / 218mm
Weight	7kg	7kg
Suitable for	PS-190 PS-191	PS-195 PS-196

06 OPERATION

REV. 2023/10

06.1. Transport and storage

Before transporting the appliance, make sure that it is correctly packaged. During transport always fix the load. If piling up is necessary, always follow the scheme shown on figure 1. Transport the load levelled and following all precepts and rules for the load transportation, ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

To unpack the device place the box on a flat, stable surface and carefully open the seal. Remove the device from the inside of the box following the scheme shown on figure 2.

Keep the original packaging for use in case of any other transport and for storage.

Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client.

The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agent which could cause damage and reduce safety characteristics.


The storage temperature must be between -20°C and 60°C.

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SIZE: 1160 x 550 x 305 mm

WEIGHT: 35 Kg.

The product must not be manipulated in any way (modification, adjustment, addition, substitution). In such cases, all responsibility for any malfunction or injury caused by the product itself will be disclaimed; furthermore, CE certification and product warranty will be considered void.

06.2. Preparation

When the product is received:

- Verify that the packaging is intact and has no signs of impacts, falls and isn't wet.
- Remove the packaging and expose the material so that all components are visible.
- Check that all the components/pieces on the section 5 are present.

The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage. In particular, check:

- General functionality of the device.
- Cleaning the device (remember that the lack of cleaning can cause the spread of infections).
- Absence of cuts, holes, breaks in the structure, including the belts.
- Correct fixing of all nuts, bolts and screws.
- Correct attachment of the belts to the chair.
- Correct operation of the belt closure.
- Condition of moving parts, wheels and belts.
- Integrity of the seams.
- There are no tubes or metal sheets that show bends or cracks.
- The backrest and the seat do not show damage or structural cracks.
- The welds are intact, without cracks or breaks.
- Wheels are securely fixed, stable and work properly.
- Wheels are free of dirt or debris.
- Brakes work properly.
- Track belts run and have the correct tension for use.
- Correct operation of the springs.
- Armrests raise and lower correctly.
- The transport handles open, close and lock correctly.
- Telescopic handles open, close and lock properly.
- Presence of all labeling.
- Lubrication of moving parts.
- The emergency vehicle is equipped with the locking system for the PS-190 chair.

If the above conditions are met, the device may be considered ready for use; otherwise it must immediately be removed from service and contact the manufacturer.

06.3. *Functioning*

ARMRESTS

Lower the armrest into the operating position, holding it and pushing it down, until you meet resistance. To close it again, lift it up to the previous position.

TRANSPORT HANDLES

To close the transport handles, press the red release button and push them down until they stop. To reopen them, lift them up until they stop.

TELESCOPIC HANDLES

To extend the telescopic grips, stand in front of the device, press the red release buttons and pull both grips to the forward position. Release the button and slide them forward slightly until they snap into the locked position.

To fold them, stand in front of the device, press the red release buttons and push both grips to the rear position. Release the button and slide them back slightly until they snap into the locked position.

BELTS

Use the belts to help support the patient in the chair.

On the backrest of the chair you will find a belt with buckle and another with anchor, and also on the seat.

You must insert the backrest belt buckle to the seat belt anchor and vice versa, until a “click” is heard, so that the belts are fastened to the patient’s torso in the shape of a cross.

Adjust the belt measurement according to the patient.

Check that the belts are properly fastened by pulling on them.

Before each use, make sure that the belts are securely attached to their respective anchor points.

FOOTREST

To open the footrest, lower it until it stops, to close it again raise it until it stops.

When using the footrest make sure that it does not interfere with the feet of the patient or the operator.

Before transferring the patient to the chair make sure the footrest is closed, open it after the patient has sat down.

When the patient has to get out of the chair, close the footrest before unfastening the belts, in order to avoid standing on it.

Keep the footrest elevated when not in use.

WHEELS BRAKE

The 2 rear wheels of the chair are equipped with brakes to prevent the chair from moving during the transfer of the patient or during a stop.

To activate them, press down on the upper part of the pedal. To remove, press down on the front end of the pedal.

Never leave the chair unattended with the patient, always keep it under control. Brakes should not be used as a substitute for operator control.

FOLDING AND UNFOLDING THE CHAIR

Before folding the chair, it is recommended to lock the rear wheels and fold the footrest, the carrying handles and the armrests, as indicated in paragraphs 03.3.1, 03.3.2, 03.3.5 and 03.3.6.

To fold the chair, locate the black release bar at the bottom of the seat. Pull the bar to the front, sliding it through the slot, until the chair is folded down. Verify that the system has been locked correctly by trying to open the chair without activating the bar, if it is locked the chair will not open.

06 OPERATION

REV. 2023/10

To unfold the chair, pull the release bar and move the backrest away from the seat until the chair is open.

Verify that the system has locked correctly by trying to close the chair without activating the bar, if it is locked the chair will not close.

Once unfolded, you can unlock the rear wheels again.

ELECTRIC TRACK

Use the electric track to facilitate the transport of the patient up stairs.

OPENING THE TRACK:

To open the track, stand behind the chair, locate the side latch, open it, and pull the black tube of the track toward you until it stops. Verify that the system has been locked correctly by pushing the track as if you wanted to close it, if it is locked the track will not close.

CLOSING THE TRACK:

To close it, stand behind the chair. Hold the two carry handles and place one foot on the red tube of the track. Press the tube down along the slot until it stops. The side latch will lock automatically. Verify that the system has been locked correctly by pulling the track as if you wanted to open it, without opening the side latch, if it is locked the track will not open.

ADJUST TRACK TIGHTNESS

Periodically check the tension of the track belts. If it is detected that the tension is not correct, it must be adjusted through the screw at the top of each strap, with the help of an Allen key.

It is important not to over-tighten the screw as it could cause the gear bearing to make noise and may even cause a failure of the mechanism.

INSTRUCTIONS FOR USE OF THE TRACK:

The chair can be operated on stairs by a single operator, but in extremely heavy patient loading procedures, in rough terrain operations and, in particular situations, the use of two operators is recommended.

Never lubricate the belts. Lubrication can cause belts to perform unpredictably, which can cause injury to the patient and/or operators.

Moisture, water, snow, ice or debris on the crawler belts or in the path of the stairs can cause uneven operation of the product causing sudden changes in the weight to be supported by the operators or destabilizing the chair.

Make sure the track and track belts are clean and dry before using the set on stairs.

Before using the track, adjust the upper handle to the desired height and fold the rear levers. Press the power button on the control panel, approach the chair to the start of the stairs and gently recline it until the tracked powertrain is on top of the first few steps.

The second operator must extend the front telescopic handles and hold them throughout the operation.

Press the UP button if you want to go up the stairs, or the DOWN button if you want to go down.

06 OPERATION

REV. 2023/10

Once the flight of stairs is finished, all the elements that are not necessary for the transfer (telescopic handles and track) must be folded back.

To stop on the stairs, you must first confirm that the track is in contact with at least three steps. You can then proceed to release the up or down button. The chair will remain stable in this position.

If your chair has a control panel that allows you to regulate the speed, you must follow this process:

- To adjust the upload speed you must press the upload button and the + or - buttons.
- To adjust the descent speed press the down button and the + or - buttons.

PANEL CONTROL

- 1 UP
- 2 DOWN
- 3 LED LIGHT
- 4 SPEED UP
- 5 SPEED DOWN
- 6 TURN ON / OFF

FIGURA 17 - FUNCIONES DEL PANEL DE CONTROL

BATTERY

Locate the slot to position the battery under the seat in the rear part. Insert the battery from top to bottom.

Once attached, you can lock it by key.

The battery charging time is 5-6 hours. When the battery is charging, the charger light is red. When fully charged, the light turns green. It is recommended to continue charging for an additional hour after the green light has come on.

The battery can last 1.5 hours once fully charged, which is equivalent to going up and down a 5-story building 30 times.

When the chair is not in use, close the battery switch.

If the chair is not used for a long time, the battery should be charged every 3 months.

The battery can be charged 500 times.

The service life of the battery is 3-5 years.

07 GENERAL MAINTENANCE

REV. 2023/10

07.1. Cleaning

It is essential to keep the equipment clean to ensure proper use and durability of the set. Thorough cleaning should be carried out periodically, especially in areas exposed to dirt that can accumulate dirt such as corners or holes.

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Do not use high pressure cleaning systems or bleach to disinfect the product, they can damage it. Instead, use a hydroalcoholic-based disinfectant, wash it with water and let it dry naturally, do not use direct heat sources to dry.

Non-compliance with cleaning operations can lead to the risk of cross infections due to the presence of secretions and/or residues. During all control and disinfection operations, the operator must wear appropriate personal protective equipment, such as gloves, glasses, etc.

When the cleaning procedure includes a disinfectant product, carefully follow the instructions of the manufacturer of the product used regarding the method of application and contact time. Make sure that you have taken all the necessary precautions to ensure that there are no risks of cross-infection or contamination of patients and operators.

MAKE SURE TO KEEP AREAS FROM FIGURES ON THE RIGHT FREE OF WATER AND HUMIDITY. TAKE SPECIAL CARE NOT TO GET THEM WET WHILE WASHING THE SET.

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals. The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

07.2. Lubricating

In general, all moving parts must be greased. Our products leave the factory completely greased and lubricated, however, it is possible that over time and use of the product the elements become degreased, either due to loss of lubrication or due to dirt.

Periodically clean and oil affected areas according to manufacturer's specifications. Check for loose, missing or worn parts.

Never lubricate the belts. Lubrication can cause the belts to perform unpredictably, which can cause injury to the patient and/or operators.

07.3. Wear areas

Periodic inspection of the different components of the system for signs of wear is a preventative measure that can reduce breakdowns. Check possible lubricant leakages, grooves or bearing in poor condition. Check also, due to intensive usage, the moving parts such the handles, the wheels and the folding system of the chair and the track.

In the event that a worn component needs to be changed, you should contact our sales department for more information on ordering spare parts and their installation.

The repair or replacement of product components must necessarily be carried out by qualified personnel, that is, who have received training from the Manufacturer, Promeba S.L. and original spare parts should always be used.

Promeba S.L. declines any responsibility for any damage, direct or indirect, resulting from improper use of spare parts and/or any repair intervention carried out by unauthorized persons.

07.4. Maintenance program

Although the frequency of the controls is determined by factors such as legal requirements, the type of use, the frequency of use, the environmental conditions during use and storage, it is placed a recommendation on the maintenance that the equipment must have for correct operation throughout its useful life.

All maintenance and revision activities must be recorded and documented with the corresponding technical intervention reports, and the documentation must be kept for at least 10 years from the end of the product's useful life and must be made available to the competent Authorities and/or from the Manufacturer, when requested.

Maintenance review	Every use	When necessary	Every month	Every year	Every 4 years
Disinfect	X				
Clean		X			
Inspect	X	X			
Grease necessary components		X		X	
Change springs				X	
Change wheels					X

For both the disinfection and cleaning procedures, point 7.1 of this manual must be followed.

During the inspection procedure, which must be done after each use, the following tasks must be carried out:

- Check that all components are present.
- Check the integrity of the device, that is, that there are no breaks, cracks, holes or cuts.
- Check the wear of the product following section 7.3 of this user guide.
- Check that the moving parts slide correctly.

07 GENERAL MAINTENANCE

REV. 2023/10

- Check that the accessories provided are present, intact and that they function correctly.

For some conditions of use, due to vibrations or impacts, certain elements may lose their tightening torque or fixing properties. Periodically check for loose items, especially moving parts. Inspect any mechanical connections. Always observe and adhere to the recommended tightening torques.

08 SPARE PARTS

REV. 2023/10

To request a spare part, indicate:

- Number and description of the set in which the necessary part is found (see table)
- Point out the piece in the drawing, or the area where it is located.

Following these instructions, our sales department will provide you with a detailed drawing of the assembly so that you can specify the reference of the required spare part.

PS1900-00000 / PS1910-00000

N°	COMPONENTS DESCRIPTION
1.01	BELT TRACK
1.02	PS-190 BATTERY
1.03	SWIVEL WHEEL PS-190 ø75mm
1.04	FIX WHEEL PS-190 ø120mm
1.05	EXTERNAL CHARGER FOR PS-190 CHAIR
1.06	ARMREST
1.07	LEFT REAR HANDLE

N°	COMPONENTS DESCRIPTION
1.08	RIGHT REAR HANDLE
1.09	REAR BRAKE FORK
1.10	FEET BELTS PS-190
1.11	SEAT BELTS PS-190
1.12	RIGHT TELESCOPIC HANDLE
1.13	LEFT TELESCOPIC HANDLE

08 SPARE PARTS

REV. 2023/10

PS1950-00000 / PS1960-00000

N°	COMPONENTS DESCRIPTION
1.01	BELT TRACK
1.02	PS-190 BATTERY
1.03	SWIVEL WHEEL PS-190 ø75mm
1.04	FIX WHEEL PS-190 ø120mm
1.05	EXTERNAL CHARGER FOR PS-190 CHAIR
1.06	ARMREST
1.07	LEFT REAR HANDLE

N°	COMPONENTS DESCRIPTION
1.08	RIGHT REAR HANDLE
1.09	REAR BRAKE FORK
1.10	FEET BELTS PS-190
1.11	SEAT BELTS PS-190
1.12	RIGHT TELESCOPIC HANDLE
1.13	LEFT TELESCOPIC HANDLE

11 LEGAL NOTICES

REV. 2023/10

- This document may contain technical inaccuracies or typographical mistakes.
- Changes are periodically added to the information herein; these changes will be incorporated in new editions of the publication.
- Promeba, S.L. reserves the right to make any modification or improvement in the products described in this publication if it is appropriate.
- Promeba, S.L. may have patents or patent applications which address themes described in this document. The possession of this document does not entitle no license to these patents. The information contained in this document does not affect or change the specifications or product specifications or warranties of Promeba, S.L.
- No part of this document shall operate as express or implied license or compensation under the intellectual property rights of Promeba, S.L. or .third party.
- All the information contained in this document has been obtained in specific environments and it is presented as an illustration. The results obtained in other operating environments may vary.
- Promeba, S.L. can use or distribute the information that the customer provides in whatever way they see fit, without incurring any obligation with the client.

12 PRODUCT WARRANTY

REV. 2023/10

Promeba, S.L. guarantees that its products have satisfactorily passed all established quality controls, both functional and material. The duration of the guarantee is 2 years from the date of purchase of the product.

This guarantee will be valid only when it is presented with the original invoice or proof of purchase (indicating the date of purchase, model and the name of the distributor) together with the defective product during the period covered by the guarantee. Promeba, S.L. reserves the right not to offer the free warranty service if the indicated documents are not presented or if the information they contain is incomplete or illegible.

1. This warranty will not apply if the model name or serial number of the product has been altered, erased, disappeared or becomes illegible.

2. This guarantee does not cover transportation costs or risks derived from the transportation of your product to and from Promeba, S.L.

3. This warranty does not cover any of the following:

- a) Periodic maintenance and repair or replacement of parts derived from normal wear and tear.
- b) Expendable material (components that are expected to need periodic replacement during the life of the product, such as belts, straps, etc.).
- c) Damages or defects derived from the improper use, operation or treatment of the product and not due to normal use of the product.
- d) Damages derived from:
 - i. Misuse, including:
 - Treatment that results in damages or physical, superficial or appearance changes of the product.
 - Installation, use or storage of the product in a way that does not

respect the instructions described by Promeba, S.L.

- Maintenance of the product in a way that does not respect the instructions of Promeba, S.L. for proper maintenance.

- Installation or use of the product in a way that does not respect the technical or safety regulations of the country where it is used or installed.

ii. Use of components not provided with the product or incorrect installation of accessory parts not previously tested.

iii. States or defects of the system in which the product is used or incorporated with the exception of other products Promeba, S.L. designed for use with the product.

iv. Use of the product with accessories, peripheral units and other products of a type, condition or standards not established by Promeba, S.L.

v. The manufacturer or distributor will be solely responsible for deciding whether to send the parts for repair, or to replace the product in its entirety. In no case will operators be sent to repair or replacement of the product.

Except in the cases mentioned above, Promeba, S.L. will make no warranties in relation to product, performance, accuracy, reliability, or fitness for logic purpose of the equipment or other type. If this exception is not legal or contemplated by current law, Promeba, S.L. will limit or exclude your warranties only to the extent permitted by applicable law.

The only obligation on the part of Promeba, S.L. in connection with this warranty is to repair or replace parts subject to the terms and conditions of this warranty.

Promeba, S.L. is not responsible for loss or damage to products, this warranty or others, including economic loss or non-assessable

12 PRODUCT WARRANTY

REV. 2023/10

damage; the price paid for the product; loss of profits, income, information, usufruct or use of the product or associated products or indirect, accidental or critical loss or damage.

This clause refers to whether the loss or damage is due to deterioration or inoperability of the associated product due to defects or unavailability of Promeba, S.L., which has caused downtime, loss of user time or business interruption.

In cases where the law prohibits or limits these liability exclusions, Promeba, S.L. will exclude or limit his liability only to the extent permitted by applicable law.

For example, there are countries that prohibit the exclusion or limitation of damages caused by negligence, reckless negligence, willful misconduct, fraud and similar actions. The responsibility of Promeba, S.L. In this guarantee will not exceed, in any case, the price paid for the product, but if the current law allows only limitations of greater responsibilities, these will apply.

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*All rights reserved. Variations can be done whitout notice.
Promeba, S.L. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.*

SCAN QR CODE

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Ctra C-16 Km 59.5 · 08650 Sallent (Barcelona) · SPAIN

T. 93 837 12 00

promeba@promeba.com · www.promeba.com