DOMUS2

Alternating Pressure Redistribution System

Operation Manual

IMPORTANT SAFEGUARDS READ ALL INSTRUCTIONS BEFORE USING

DANGER - To reduce the risk of electrocution:

- 1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store this product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING - To reduce the risk of burns, electrocution, fire or injury to persons:

- 1. Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
2. This system is not for use with patients who have a spinal cord injury.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.

- 4. Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Apex Medical Corp. for examination and repair.

- 6. Keep the cord away from heated surfaces.
7. Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.

- 8. Never drop or insert any object into any opening or hose.
9. Do not modify this equipment without authorization of the manufacturer.
10. Mattress covers have passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.

- 11. Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.

CAUTION - If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

NOTE, CAUTION AND WARNING STATEMENTS:

- NOTE - Indicate some tips.
CAUTION - Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property
WARNING - Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.


1. INTRODUCTION

This manual should be used for initial set up of the system and for reference purposes.

1.1 General Information

The system is a high quality and affordable mattress system suitable for treatment and prevention of pressure ulcers.

The system has been tested and successfully approved to the following standards:

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2:2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
Increase the separation between the equipment.
Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
Consult the manufacturer or field service technician for help.

1.2 Intended Use

This product is intended:

- to help and reduce the incidence of pressure ulcers while optimizing patient comfort.
for long term home care of patients suffering from pressure ulcers.
for pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure ulcers.

NOTE: This equipment is not suitable for use in the presence of flammable anesthetic mixture with air or with pure oxygen or nitrous oxide.

1.3 Warranty-

The Company warrants the pump at the time of its original purchase and for the subsequence time of period of one year.

The company warrants the bubble pad at the time of its original purchase and for the subsequence time of period of six months.

The warranty does not cover the followings:


- 1. The series number label of pump or pads is turn off or cannot be recognized.
2. Damage to the pump or bubble pad resulting from misconnections with other devices.
3. Damage to the device resulting from accidents.

SYMBOL DEFINITION

Manufacturer, Authorized representative in the European community, Attention, should read the instructions, Catalogue Number, Class II Equipment, "BF" symbol, indicate this product is according to the degree of protecting against electric shock for type BF equipment, IP21, Waste Electrical & Electronic Equipment (WEEE): This product should be handed over to an applicable collection point for the recycling of electrical and electronic equipment.

2. PRODUCT DESCRIPTION

2.1 Pump Unit

2.2 Front panel

- 1. Pressure Adjust Knob
2. Main Power Switch

3. INSTALLATION

Unpack the box and check the package contents for completeness. If there are any damages, please immediately contact your supplier or Apex Medical Corp

3.1 Pump & Mattress Installation

- 1. Place the DOMUS 2 on top of a foam or mattress bed. Secure the mattress firmly by fixing the straps to the bed frame if available.
2. Hang the pump onto the footboard and adjust hangers so the pump is secured in an upright position; or place the pump on a flat surface.
3. Connect air hose connectors from air mattress to the pump unit.
4. Plug the power cord into electrical outlet.

CAUTION: The pump should only be used with mattress recommended by the manufacturer. Do not use it for any other purpose.

- Turn the main power switch found from the right side of the pump to ON position.

Several installation tips are listed below:

After installation, the extra length of the power cord, if any, should be neatly arranged to avoid any tripping accidents. The EQUIPMENT should be firmly placed at position where users/doctors can access easily.

4. OPERATION

NOTE: Always read the operating instruction before use.

4.1 General operation

- Switch on the main power switch on right side of the pump.
Adjust the pressure setting based upon patient comfort level by turning the pressure adjustment knob clockwise to increase firmness.
NOTE: Every time the mattress is set up for use, it is recommended that the pressure first to be set to the max. The user / career can then adjust air mattress weight levels to the desired softness after set up has been completed

4.2 Emergency operation

When there is a need to perform emergency CPR on the patient, pull the quick release CPR tag located at the head-end of the mattress on the right hand side if available. The quick connector found from the pump unit can also be disconnected to hasten deflation. Be sure to reconnect the quick connector to the pump unit once restore the power supply.

5. CLEANING

It is important to follow the cleaning procedures before using the equipment on human bodies; otherwise, patients and/or doctors may have the possibility of getting infection.

CAUTION- Do not immerse or soak pump unit.

Wipe the pump unit with a damp cloth and a mild detergent. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit.

Wipe down the mattress with warm water containing a mild detergent. The cover may also be cleaned by using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.

CAUTION - Do not use phenolic based product for cleaning.
CAUTION - After cleaning, dry the mattress without direct exposure sunlight.

6. STORAGE

- 1. Lay the bubble pad or mattress on a flat surface and upside down.
2. Roll-up the mattress from the head end towards the foot end.
3. Foot-end strap can then be stretched around the rolled pad/ mattress to prevent unrolling.

NOTE: Do not fold, crease or stack the mattresses.

7. MAINTENANCE

- 1. Check main power cord and do not plug it if there is an abrasion or excessive wear.
2. Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
3. Check the airflow from the air hose. The airflow should alternate between each connector every half-cycle time.
4. Check the air hoses if there is kink or breaks. For replacement, please contact Apex Medical Corp. or your suppliers.

8. EXPECTED SERVICE LIFE

The products are intended to offer safe and reliable operation when used or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function. Otherwise, service and inspection of the devices generally should not be required.

9. TROUBLESHOOTING

Table with 2 columns: Problem, Solution. Rows include Power is not ON, Patient is bottoming out, Mattress form is loose, No air produced from some air outlets of the air tube connector.

NOTE: If the pressure level is consistently low, check for any leakage (tubes or air hoses). If necessary, replace any damaged tubes or hoses or contact your local qualified dealer for repair.

10. TECHNICAL SPECIFICATION

Technical specification table for Domus 2 (9P-047580) and Domus 2 (9P-047500). Includes Model, Power Supply, Fuse Rating, Dimension, Cycle Time, Weight, Environment (Atmospheric Pressure, Temp, Humidity), Classification, and Mattress Specification.

NOTE: Please follow national requirements to dispose the unit properly.

Appendix A: EMC Information

Guidance and Manufacturer's Declaration- Electromagnetic Emissions: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Table with 3 columns: Emissions Test, Compliance, Electromagnetic Environment-Guidance. Rows include RF emissions CISPR 11, Harmonic emissions IEC61000-3-2, Voltage fluctuations / Flicker emissions IEC61000-3-3.

Guidance and Manufacturer's Declaration- Electromagnetic Immunity: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Table with 4 columns: Immunity Test, IEC60601 test level, Compliance, Electromagnetic Environment-Guidance. Rows include Electrostatic Discharge (ESD), Electrical fast transient/ burst, Surge, Voltage dips, short interruptions and voltage variations on power supply input lines, Power frequency magnetic field.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity: This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Table with 4 columns: Immunity Test, IEC60601 test level, Compliance, Electromagnetic Environment-Guidance. Rows include Conducted RF, Radiated RF.

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

Recommended separation distances between portable and mobile RF communications equipment and this device: This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Table showing separation distance according to frequency of transmitter. Columns: 150 kHz to 80 MHz, 80 MHz to 800 MHz, 800 MHz to 2.5 GHz. Rows: 0.01, 0.1, 1, 10, 100 W.

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Logos for APEX MEDICAL S.L., APEX MEDICAL CORP., and Apex Medical (Kunshan) Corp. with addresses and website information.

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