

LUXURY SWIVEL SEAT LEATHER

EASY ACCESS EVERY TIME



movinglife[®]

SEATING & SUPPORT



ATTO /ATTO
SPORT Armrests



Seat cushion



Seat cushion
leather



Luxury Swivel
Seat Leather



Seat & Cushion
Support Belt



Luxury Orthopedic
Backrest

STORAGE



Essentials
Pouch



Front & rear
basket



Caryall
& Cushion



Caryall
& Cushion leather



Pet Carrier &
Cushion Black



Flight Kit



Oxygen Bottle
Holder & Case

TECHNICAL



Cane / Crutches
Holders




Mobile Phone
Holder



Rear Light Lifting
Handle



ATTO Left Hand
Throttle Kit



Standard Flight
Battery



Launch - H2 2023

Long distance
battery - not flight
approved



Cable kit


ATTO Accessories are created specifically to enhance, support and work with the ATTO mobility scooters in all user modes; while driving, while wheeling along as trolley and when folded for storage or transportation.

SCAN TO LEARN
HOW TO COMBINE
OUR ACCESSORIES




How to Use

LUXURY SWIVEL SEAT LEATHER




01

Fit the cushion over the bottom of the seat and pass the connecting strap through the gap in the seat



02

Fasten the connection strap to the bottom of the cushion



03

The product is ready to use. Must be removed before folding



04

Attach to ATTO Carryall & Cushion dedicated upper zip