Technical information Illico



* weight excl. side rails, lifting pole and triangle. Deviations in sizes and weight (15 mm and 1,5 kg) are possible

Different boards, different side rails



To match your interior, the **Illico** bed can be equipped with oxide boards or beech boards.



The metal siderails, one per side are easy to mount and to lower with just one pull. The wooden side rails, 2 per side, are complete with plastic covers at the ends and are easy to mount.

Your guarantee of quality

Regarding quality, the Illico bed complies with the highest norms laid down by certified test institutes in most European countries.

N.V. Vermeiren N.V.

Vermeirenplein 1-15 B-2920 Kalmthout Belgium

Tel.: +32 (0) 3 620 20 20 Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

VERMEIREN

Opening up new horizons in mobility



The bedframe is also available in **ABS** material for easy disinfection.



For more information, contact your local dealer



Images for illustration purposes only. Certain features may be optional

Opening up new horizons in mobility

Subject to technical modifications.

Illico Electrical homecare bed

Comfortable and practical

- Electrical height adjustment mattress support 298 mm to 828 mm
- Electrical adjustment with the included controller
- Electrical head adjustment
- Electrical knee buckle, the foot part can be adjusted manually
- Mattress size 900 mm x 2000 mm
- Handset with lock function
- Maximum load 170 kg (patient and mattress)
- Two pair of siderails
- Selfraising device can be placed left or right
- All 4 wheels can be blocked separately.
- Different front and rear wooden boards
- Easy to (dis)assemble
- Deliverable transportset



Straightforward handset

with locking functions



Easy removal of the boards and the siderails





wheels with locks

Big, star shaped, screws for easy (dis)assembly of the bed frame

The lifting pole can be positioned on the left or the right side of the bed. The strap can be adjusted to the desired length.

4 mattress positioners (2 at the head part and two at the foot part) keep the mattress in perfect place.

A single press on the button releases the side rails.

all set.

The handset comes with a long extension cable so it can be placed anywhere on the bed.

The foot part can be raised manually into different positions.



A transportkit on wheels for smooth storage



Disconnecting the bed frame without additional tools



Ultra EASY assembly without extra tools

No more searching for screwdrivers or Allen keys with our **Illico** bed. It takes just one person to assemble or disassemble this bed in a matter of minutes. Remove the side rails, the wooden boards, disconnect the bed frame and you are







